



此模具拥有切割塑料、纸张和橡胶样品，制成样品后，做拉伸、撕裂试验等。

应用:

- 瓦楞板
- 泡沫
- 纸张
- 纸张
- 塑料

产品:

- 可反复再磨快

参考标准:

- BS 2782-3
- AS 1683.11
- AS 2282.7
- ASTM D470
- ASTM D1004
- ASTM D1822
- ASTM D638
- ISO 37
- ISO 6259
- K6301
- AS 1301.404
- AS 1180.2
- ASTM D412
- ASTM D624
- ASTM 1660.2
- ASTM D3574
- AS 1145
- ISO 527
- K6251
- AS/NZS 1660.2.3



可选压头

可选择配套使用产品



手动切压机




气动切压机




液压切压机

All dimensions in mm




C0024 – Type A

- ISO 37-Type 2
- ISO 527-2-Type 5A
- AS 1145-Type 5A
- K6251-Type 6
- AS/NZS 1660.2.3




C0024 – Type B

- ASTM D412-C
- ISO 6259-Type 2
- K6251-Type 5
- ISO 37 – Type 1
- ISO 527-3 – Type 5
- BS 2782-3 (Method 320A)



C0024 – Type C


- ASTM D412-D



C0024 – Type D


- JIS K6301-Type 3
- K6251- Type 3
- ISO 37 – Type 1A

All dimensions in mm




C0024 - Type E

- ASTM D638-02a-Type 4
- ASTM D6693




C0024 - Type F

- JIS K6301-Type 1



C0024 - Type G


- JIS K6301-Type 2



C0024 - Type H


- JIS K6301-Type 4

All dimensions in mm




C0024 – Type I

- ISO 37 – Type 3
- AS/NZS 1660.2.3
- DIN 53504




C0024 – Type J

- ISO 37 – Type 4
- ISO 527-2 – Type 5B



C0024 – Type K


- K6301- B
- ASTM D624 – Type C
- ASTM D1004-03



C0024 – Type L


- ASTM D638-02a-Type 1

All dimensions in mm




C0024 – Type M

- ASTM D638-02a-Type 2




C0024 – Type N

- ASTM D638-02a-Type 3



C0024 – Type O


- ASTM D638-02a-Type 5



C0024 – Type P


- DIN ISO 34-1
- K6301-Type A
- ASTM D624 – Type B

All dimensions in mm




C0024 – Type Q

- ASTM D1822-Type L




C0024 – Type R

- ASTM D1822-Type S



C0024 – Type S


- ASTM D3574 – Type E



C0024 – Type T


- ASTM D3574 – Type F

All dimensions in mm




C0024 – Type U

- ISO 527-2 – Type 1A




C0024 – Type V

- ISO 527-2 – Type 1B
- ISO 527-3 – Type 1B
- BS 2782-3 (Method 320B)
- BS 2782-10 – Type 1




C0024 – Type W

- ISO 527-2 – Type 1BA



C0024 – Type X


- ISO 527-2 – Type 1BB



C0024 – Type Y


- ISO 527-3 – Type 2

All dimensions in mm




C0024 - Type Z

- ISO 527-3 - Type 4




C0024 - Type Aa

- ASTM D412 - Type A



C0024 - Type Ab


- ASTM D412 - Type B



C0024 - Type Ac


- ASTM D412 - Type E

All dimensions in mm



C0024 – Type Ad

- ASTM D412 – Type F



C0024 – Type Ae

- ISO 1798

* Other custom dies are available upon request

All dimensions in mm