

TESTING INSTRUMENTS FOR
A MEASURABLE DIFFERENCE...

IDM
instruments

CUTTING DIES – RULER TYPE

MODEL: IDM-C0025-M1



INTRODUCTION

The Ruler Type Cutting Dies are used for cutting of plastic film, paper and rubber samples used for tensile and tear testing (dumbbell, trouser, etc). They are manufactured with a ply backed knife-edge, and can be manufactured to suit various international test standards.

APPLICATIONS

- Plastic Film
- Paper
- Rubber

FEATURES

- Low volume work

STANDARDS

- BS2783
- AS1301.404
- AS1683.11
- AS1180.2
- AS2282.7
- ASTM D412
- ASTM D470

- ASTM D624
 - ASTM D638
 - ASTM D1004
 - ASTM 1660.2
 - ASTM D1822
 - ASTM D3574
- ** Custom Sizes Available

RELATED ITEMS



Arbour Cutting Press
(IDM-C0022-3-M1)



Tool Steel Cutting Dies
(IDM-C0024-M1)



Pneumatic Cutter
(IDM-C0050-M1)

WARRANTY AND CALIBRATION SERVICES

- 1 year Warranty

Our Preventive Maintenance and Calibration (PM&C) program has been designed to make the maintenance and calibration of your valuable testing equipment more cost effective by preventing breakdowns and downtime by regular calibration, service and replacement of defective parts. Take to us about this today.