



The **PS200** polishing machine is used to polish metal material, can get a bright metal surface as a mirror, then it can be used for observe and determine the metallurgical structure under the microscope.



The **PD200** is a double-disk table polishing machine, which can operate two samples simultaneously.



Both units offer good stability and operation with low noise.

Specifications:

MODEL	PS200	PD200
Work Disc No.	1	2
Disc Diameter	200mm	
Speed	1400rpm	
Power	3P / 380V / 50Hz or 3P / 220V / 60Hz	
Dimensions	360 x 450 x 320mm	730 x 405 x 340mm
Weight	16kg	35kg
Accessories Included	1 x Polishing Cloth	2 x Polishing Cloths

Other Applicable Units: Please see MP Series Metallographic Sample Grinder/Polishers

