Grammage Cutter Model: C0042



The Grammage Cutter has been designed to cut 100cm² circular paper samples to assess Grammage.

Grammage in paper products is calculated from the ratio of mass to area. This is an important factor to know as most paper is bought and sold in accordance with its mass per unit area.

Other values such as bursting strength, thickness and bulk are all determined and interpreted using knowledge of its Grammage.

Applications:

- Paper
- Paperboard
- Solid Fibreboard

Features:

- Sample support tray
- Bench mounting
- Easy operation

Benefits:

- Easy to use
- Fast results
- Accurate

Options:

N/A

Connections:

N/A

Dimensions:

H: 247mmW: 400mm

• **D**: 410mm

• Weight: 30kg

