

Light Transmission Tester

Model: L0005

a measurable difference...

IDM[®]
instruments

The Light Transmission Tester has been specially designed and manufactured to evaluate the light transmission properties of various shaped plastic containers. Slight variations in the mix of raw materials used for plastic containers often leads to dramatic differences in the look of the product on the shelf.

The IDM Light Transmission Tester provides a digital readout allowing you to quantify samples based on your in-house pass/fail criteria, and obtain consistent results while minimising your inputs.

Applications:

- Plastic containers

Features:

- Light Meter
- 15W Light Globe
- Spring Collar to maintain pressure on sample
- Bench-top operation
- Quick and simple process

Benefits:

- Easy to use
- Fast results
- Accurate

Standards:

- N/A

Connections:

- **Electrical:** 220/240 VAC @ 50 HZ or
110 VAC @ 60 HZ
(please specify when ordering)

Dimensions:

- **H:** 455 mm
- **W:** 143 mm
- **D:** 230 mm
- **Weight:** 9 kg

