



The Random Tumble Pilling Tester performs tests on fabric under free roll friction conditions.



### Standards:

- ASTM D1375
- ASTM D3512
- GB/T4802.4
- ISO 12945.3

### Applications:

- Textiles
- Woven Fabrics

### Specifications:

<b>No. of Test Areas (Boxes)</b>	4
<b>Drum</b>	(146±1) x (152±1) mm
<b>Cork Lining (LxWxH)</b>	(152±1) x (146±1) x (1.5±1) mm
<b>Stirring Rod Length</b>	121 mm
<b>Rotating Speed</b>	1200±1 r/min
<b>Noise</b>	≤50db
<b>Compressed Air</b>	0.014-0.021 Mpa
<b>Power</b>	AC220V, 50Hz, 80W
<b>Dimensions</b>	480 x 400 x 680 mm