

The Automatic Notch Machine is mainly used to notch on specimens for IZOD impact test and CHARPY impact test.

ZQK-20 automatic notch machine can prepare Izod and Charpy impact test samples (types A, B, and C).

Other types of samples can be notched through changing milling cutter and standard specimen thickness profile.

Applications:

- Izod Test Specimens
- Charpy Test Specimens

Features:

- Impact Angle: 90°
- Capacity: 0.25/0.5kg/cm
- Resolution: 0.005kg/cm

Standards:

- GB/T1043-1993
- GB/T1843-1996
- ISO 179-1993
- ISO 180-1993

Connections:

- Electrical: 220V, 50/60Hz

Dimensions:

- W: 340 mm
- L: 700 mm
- H: 600 mm
- Weight: 91 kg

Benefits:

- Easy to use
- Accurate

