

Cavity Sample Moulds

Model: D0004 & D0005

a measurable difference...

IDM[®]
instruments

These cavity sample moulds, available in both 2 and 4 cavity sizes, are constructed from hardened stainless steel and allow for plastic and rubber samples to be made with ease, repeatability and confidence.



Applications:

- Plastic
- Rubber
- Synthetic Rubber

Features:

- 2 or 4 Cavities
- Hardened stainless steel
- **For 400mm x 400mm platen:**
150mm x 150mm x 2mm deep cavities
- **For 300mm x 300mm platen:**
110mm x 110mm x 2mm deep cavities

Standards:

- ASTM D3182

Options:

- Custom made cavity sizes
(Available on request)

Connections:

- N/A

Dimensions:

- **H:** 20mm
- **W:** 380mm
- **D:** 380mm
- **Weight:** 10kg